

Pair the Amwell app to your Tyto Device

For your Tyto Device to communicate with the Amwell app, you need to pair it. You will use your home Wi-Fi or other Wi-Fi whenever using the Tyto device with Amwell, so please make sure the smartphone or tablet you are using is connected to your home Wi-Fi and that you are in a location where the signal is strong (e.g., close to the router).

Please have your home Wi-Fi password available as you will need it to complete the pairing process.

It is always preferred that you use your home Wi-Fi Internet connection when using the Tyto device, but if you need a backup or if you are traveling, you may use your smartphone as a hotspot Internet connection.


Step 1

Download the Amwell Doctor's Visit 24/7




Step 2

Log into the Amwell app.
Set up account if haven't already
Amwell will try to set up account if we have an email to share with them.




Step 3

Select Mohawk Pediatric Practice




Step 4

Select physician.
List would be pediatricians only for this practice.




Step 5

Check the reason for visit and start pairing process. Please make sure the Tyto device is turned on.



Step 6


Enter name of Wi-Fi network and password (once connected, you will only change if on a different network or hot spot (if applicable). Hit continue.



Step 7

Use the camera on the Tyto Device to scan the barcode on the Amwell as shown.

Next. If this screen is not displayed automatically, press the red pairing icon, located in the upper right corner of the App and the Device pairing screen will be displayed.



Step 8

Device will connect to Wi-Fi and show correct pairing once it is complete.

Complete intake instructions.
Cost is \$20.